



ENDEAVOR™ 5000 CROSS REFERENCE

Light Duty Orbital Shakers



e-E51LD0403



e-E51LD0420



e-E51LDMP03

SHLD Light Duty Shakers to Endeavor 5000

Old SHLD Models				NEW Endeavor 5000			
Model	Control	Orbit	Capacity	Model	Control	Orbit	Capacity
SHLD0415AL	Analog	15 mm	3.6 kg	e-E51LD0420	Digital	20 mm	4.0 kg
SHLD0415DG	Digital	15 mm	3.6 kg				
SHLD0403Dg	Digital	3mm	3.6 kg	e-E51LD0403	Digital	3 mm	4.0 kg
SHLDMP03DG	Digital	3 mm	4 × Microplates	e-E51LDMP03	Digital	3 mm	4 × Microplates, 4 × Deepwell Plates

Endeavor 5000 Feature Enhancements

- Increased capacity with the ability to shake up to 6 × 250 mL flasks or 4 × 500 mL flasks on the standard tray
- Ability to use an optional universal platform to expand capacity up to 3 × 1 L fasks
- Programmable—save time and guarantee reproducible results
- RS232 connection allows for datalogging and remote functionality
- Increased durability with IP21 Rating

SHLD Light Duty Shakers			NEW Endeavor 5000 Light Duty Shaker		
30391891	SHLD0415AL, US	15 mm Orbit, Analog Control, 120V	30770910	e-E51LD0420 AM	20 mm Orbit, Digital Control, 120V
30391898	SHLD0403DG, US	3 mm Orbit, Digital Control, 120V	30770909	e-E51LD0403 AM	3 mm Orbit, Digital Control, 120V
30391905	SHLDMP03DG, US	Microplate 3 mm Orbit, Digital Control, 120V	30770911	e-E51LDMP03 AM	Microplate 3 mm Orbit, Digital Control, 120V
30391912	SHLD0415DG, US	15 mm Orbit, Digital Control, 120V	30770910	e-E51LD0420 AM	20 mm Orbit, Digital Control, 120V
30391893	SHLD0415AL, EU	15 mm Orbit, Analog Control, 230V	30770913	e-E51LD0420 EU	20 mm Orbit, Digital Control, 230V
30391900	SHLD0403DG, EU	3 mm Orbit, Digital Control, 230V	30770912	e-E51LD0403 EU	3 mm Orbit, Digital Control, 230V
30391907	SHLDMP03DG, EU	Microplate 3 mm Orbit, Digital Control, 230V	30770914	e-E51LDMP03 EU	Microplate 3 mm Orbit, Digital Control, 230V
30391914	SHLD0415DG, EU	15 mm Orbit, Digital Control, 230V	30770913	e-E51LD0420 EU	20 mm Orbit, Digital Control, 230V